










# AED COMPARISON CHART

| AED Model  | Dust/Water Resistance | Size (inches)     | Weight | Battery Life        | Pads Life | Child Defib Method            | CPR Support/Feedback  | Warranty |
|--|-----------------------|-------------------|--------|---------------------|-----------|-------------------------------|---|----------|
| <a href="#">Physio LP CR2</a><br>             | IP55                  | 8.0 x 7.25 x 1.9  | 2.4    | 4 years             | 4 years   | Button<br>(built in function) | Audio ✓ Metronome   | 8 years  |
| <a href="#">ZOLL AED+</a><br>                 | IP55                  | 5.25 x 9.5 x 11.5 | 6.7    | 5 years             | 5 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 7 years  |
| <a href="#">ZOLL AED3</a><br>                 | IP55                  | 5 x 9.3 x 9.7     | 5.5    | 5 years             | 5 years   | Button<br>(built in function) | Audio ✓<br>Real-time CPR Feedback<br>Visual ✓ Real-time meter | 8 years  |
| <a href="#">HeartSine SAM 350P</a><br>        | IP56                  | 8 x 7.25 x 1.9    | 2.4    | 4 years             | 4 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 8 years  |
| <a href="#">HeartSine SAM 360P</a><br>       | IP56                  | 8 x 7.25 x 1.9    | 2.4    | 4 years             | 4 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 8 years  |
| <a href="#">HeartSine SAM 450P</a><br>      | IP56                  | 8 x 7.25 x 1.9    | 2.4    | 4 years             | 4 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 8 years  |
| <a href="#">Cardiac Science G5</a><br>      | IP55                  | 3.4 x 9.0 x 11.8  | 5.7    | 4 years             | 2 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 8 years  |
| <a href="#">Defibtech Lifeline</a><br>      | IP54                  | 8.5 x 11.8 x 2.7  | 4.2    | 5 or 7 year options | 2 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome<br>Visual ✓                                 | 8 years  |
| <a href="#">Defibtech Lifeline View</a><br> | IP54                  | 7.3 x 9.5 x 2.3   | 3      | 4 years             | 2 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 8 years  |
| <a href="#">Philips FRx</a><br>             | IP55                  | 2.4 x 7.1 x 8.9   | 3.5    | 4 years             | 2 years   | Child Key<br>(added cost)     | Audio ✓ Metronome   | 8 years  |
| <a href="#">Philips OnSite</a><br>          | IP21                  | 2.8 x 7.4 x 8.3   | 3.3    | 4 years             | 2 years   | Pedi Pads<br>(added cost)     | Audio ✓ Metronome   | 8 years  |